S P *E* C S[®]

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
USA

**SPECS devices operate
with lethal high-voltage.**

**All SPECS safety instructions
are to be strictly followed.**

**SPECS safety instructions
are to be kept in close proximity
to the device, so that they will be
easily and quickly
accessible for all persons.**

Correct Operation and use only for the intended purpose will ensure that the device will give trouble-free performance and remain in safe condition.

Regulations for protection against radiation sources must also be observed.

INFORMATION:

For more detailed instructions on maintenance and for requirements governing safety of the device, the relevant national regulations are to be noted and strictly observed.

If these regulations are unknown, a specialist is to be consulted in all cases!

Surrounding conditions

- After the device is unpacked, it must be let stand for at least 24 hours at the installation location (possible formation of condensate water).
Only then should the power connection be made.
- The device must then be operated in a dry, covered area.
- The device must be placed on a flat level and no-slip worktop and be fitted into a suitable 19" rack.
- The device must not be exposed to extremely high or low temperatures.
- The device must not be exposed to direct sunlight and not be stored, installed or operated close to heating or heat generating devices and air-conditioners.
- The admissible range for operation temperature is defined between 5°C and 40°C.
- The device must not be operated in areas with high humidity (for example frost, dew, condensation or rain water), draught, dust accumulation or danger from oil splashes, iron dust and corrosive liquids.
- The maximum relative air humidity is defined at 80% at 31°C (50% at 40°C).
- The device must be placed so that the air ventilation is not obstructed. The specific instructions in the operation instructions are to be followed.
- Trouble-free operation can be warranted for a height of up to 2000 meters above sea level.
- The device corresponds to protection Class 1 with protective conductor /IEC 335-1 /EN 60335-1.

Danger for pacemakers

- Persons with cardiac pacemakers and other surgical implants should not expose themselves to electromagnetic fields unless their pacemakers or implants comply with DIN VDE 0750 section 9A, EN 50061, ISO 601 -2-31 * CEI.

Mains connection

visual
check

visual
check

- Care is to be taken that the voltage in the mains agree with the specifications on the identification plate on the rear of the device or in the fuse window / panel.
- Mains power supply-voltage fluctuations must not exceed $\pm 10\%$ from the rated value.
- Mains power supply-frequency fluctuations must not exceed $\pm 10\%$ from the rated value.
- The device must under no circumstances be connected with mains power supply outlets already used by other devices which could (possibly) cause interference.
- Before the device is switched on, all plugs must be securely inserted into the corresponding jacks. Screw connections are to be securely tightened and bayonet connections must lock in.
- If the device emit sparks, smokes or smells of burning, it must be switched **off** immediately and / with all peripheral devices.
- Any **EMERGENCY OFF** device fitted should also be activated.
- **Never use water to extinguish fire!**

Operation

- Before the device is started, the operator is to ensure that it is in the correct condition and operating safely.
- Do not remove safety fittings, such as covers, screw joints, locking wires and seals for / before operation.
- Use the device only for the correct purpose.


Switching off

- After the device has been switched off, the capacitors may continue to be charged with lethal voltage for some time. Check voltage in **all** capacitors before any maintenance work is commenced.
- Work with radiation sources, cables, connected peripheral devices or cables only when the sources and all peripheral devices have been switched off.
- After the device has been switched off, cable connections may be removed only after a waiting period of one minute has elapsed.

Opening the device


- The device may only be opened by specialists with detailed knowledge of the devices / systems on the basis of special training and experience.
- Before the device is opened, all terminals must be disconnected from the power supply:
Press the power switch (OFF or O Position) and remove the power plug (only on plug casing).
- In the case of devices with voltage over 5kV, covers are only to be opened with suitable tools. Care is to be taken that only non-conductive, high voltage tools must be used.
- **Never put a device into Operation when the casing is open.**

Maintenance

visual
check

- Devices, connections, power plugs and measuring equipment are to be subject to inspection for noticeable defects every six months.
- Cover all high-voltage conducting parts during maintenance work on opened devices.
- Acoustic and visual warning Signals in the area of the devices / units / Systems are to be heeded.
- As protection against hazardous electric shocks, external earthing equipment (twisted copper wire of 10mm to 15mm in width) must be subject to regular mechanical and electrical inspections for corrosion. The ground wire should be examined for power resistance every three months.
- If the device / system is operating in an unknown operating condition, it should only be approached with the utmost caution:
Left hand is to be slipped into the trouser pocket or behind the belt.
Right hand carefully approaches the device.

Cleaning


- The device may only be cleaned by qualified personnel. Before cleaning the device, care is to be taken that all peripheral devices are switched off and all plugs are removed from the jacks (complete disconnection of mains terminals).
- Only a soft dust-free cloth may be used. The cloth should be lightly moistured with water or a mild, non-aggressive liquid detergent (observe instructions for use).
- Specific care is to be taken that no moisture enters the inside of the device.
- **Never** use strong chemical agents, such as alcohol, paraffin or paint dilution.


Accessories

- Only use original cables, plugs and accessories manufactured by SPECS. These components are to be subject to regular tests. Care is to be taken that they show no evidence of any mechanical and / or electrical defects. In case of doubt, they are to be replaced immediately by original SPECS components.
- When fuses are being replaced, only the same type must be used: voltage strength, strength of rated current and switch characteristics.
- If cables, radiation sources (ions, x-rays and electrons) and / or other devices not manufactured by SPECS are connected and put into operation or used, there can be no warranty for the correct function for connected devices and systems nor for the safety of personnel.
- If there is any doubt, SPECS customer service is to be consulted immediately.

Environment protection

- Replaced devices, component groups or elements may not be disposed of as normal household refuse. These are special wastes and should be disposed of according to their composition so that no damage is caused to the environment. They can be returned to the manufacturer or distributor.

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All rights reserved. No part of this manual may be reproduced without the prior permission of SPECS® GmbH.